

Product Datasheet

Rab4 antibody (orb151572)

Description	Rabbit polyclonal to Rab4 (ATTO-488). Rab4 is a 25kDa member of the Rab family of small guanosine tr...
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	ATTO-488
Tested Applications	ICC, IHC, WB
Immunogen	C-terminal peptide from human Rab4
Target	Rab4
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:1000 dilution of SPC-141 was sufficient for detection of Rab4 in 10 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	26kDa
Uniprot ID	P20338
NCBI	NP_004569.2
Entrez	5867
Dilution Range	WB (1:1000), IHC (1:100), ICC/IF (1:150)
Expiration Date	12 months from date of receipt.