

Product Datasheet

CDC37 antibody (orb151597)

Description	Rabbit polyclonal to CDC37 (APC). Hsp90 co-chaperone Cdc37 is a protein that is encoded by the CDC37...
Species/Host	Rabbit
Reactivity	Human
Conjugation	APC
Tested Applications	ICC, WB
Immunogen	Native human Cdc37, full length
Target	CDC37
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1.68 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	A 1:2000 dilution of SPC-142 was sufficient for detection of cdc37 in 20 µg of HeLa cell lysate by ECL immunoblot analysis.
Clonality	Polyclonal
MW	44.5kDa
Uniprot ID	Q16543
NCBI	NP_008996.1
Entrez	11140
Dilution Range	WB (1:2000), ICC/IF (1:200)
Expiration Date	12 months from date of receipt.