

Product Datasheet

Human /mouse/Rat TGF-Beta 3 Protein [Out of stock] (orb1516319)

Description	Recombinant Human /mouse/Rat TGF-Beta 3 Protein
Endotoxins	< 1.0 EU/μg as determined by LAL
Target	TGFB3
Form/Appearance	Lyophilized powder
Concentration	> 0.5 mg/ml
Storage	Stored at -70°C or -20°C. Avoid repeated freeze/thaw cycles.
Tag	Tag free
Note	For research use only
Application notes	Centrifuge tubes before opening. Reconstituting to a concentration more than 100 μg/ml is recommended. Dissolve the lyophilized protein in 4 mM HCl.
Protein Sequence	300-412/412
Purity	> 95% as determined by Tris-Bis PAGE; > 95% as determined by SEC-HPLC
MW	Due to glycosylation, the protein migrates to 13-15 kDa based on Tris-Bis PAGE result. kDa
Uniprot ID	P10600
Expiration Date	6 months from date of receipt.