

Product Datasheet

GPC3 Antibody (orb1517320)

Description	GPC3 Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	Recombinant Human Glypican-3 protein (25-559AA)
Target	GPC3
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	P51654
Expiration Date	12 months from date of receipt.