

Product Datasheet

SLN Antibody (orb1517331)

Description	SLN Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IHC
Immunogen	Recombinant Human SLN protein (1-31AA)
Target	SLN
Preservatives	PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
Uniprot ID	O00631
Expiration Date	12 months from date of receipt.