

Product Datasheet

FILIP1L Antibody (orb1517346)

Description	FILIP1L Monoclonal Antibody
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Human Filamin A-interacting protein 1-like protein (1-200AA)
Target	FILIP1L
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q4L180
Expiration Date	12 months from date of receipt.