

Product Datasheet

CNR1 Antibody (orb1517450)

Description	CNR1 Monoclonal Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human CNR1 protein
Target	CNR1
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	hIgG4(S228P)
Clonality	Monoclonal
Uniprot ID	P21554
Expiration Date	12 months from date of receipt.