

Product Datasheet

NECTIN4 Antibody (orb1517477)

Description	NECTIN4 Monoclonal Antibody
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human NECTIN4 protein
Target	NECTIN4
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	hIgG1
Clonality	Monoclonal
Uniprot ID	Q96NY8
Expiration Date	12 months from date of receipt.