

Product Datasheet

ApoE Antibody (orb1517484)

Description	ApoE Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Agarose
Target	ApoE
Preservatives	75 mM PBS, 75 mM NaCl, 0.5 mM EDTA, 0.02% NaN ₃ , 0.1mM PMSF, pH 7.3
Form/Appearance	Whole IgG
Concentration	5 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	12 months from date of receipt.