

Product Datasheet

SARS-CoV-2 S-RBD(319-541)/HIS Antibody (orb1517534)

Description	SARS-CoV-2 S-RBD(319-541)/HIS
Conjugation	Unconjugated
Tested Applications	ELISA
Target	SARS-CoV-2 S-RBD(319-541)/HIS
Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Concentration	1000 µg/ml
Storage	Store at -20°C to -80°C
Note	For research use only
Clonality	Recombinant
Uniprot ID	P0DTC2
NCBI	YP_009724390.1
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	6 months from date of receipt.