

Product Datasheet

BTBD12 Antibody (orb1517727)

Description	BTBD12 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 1784 and 1834
Target	BTBD12
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	250 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8IY92
NCBI	NP_115820.2
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.