

Product Datasheet

G9A/EHMT2 Antibody (orb1517994)

Description	G9A/EHMT2 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 1160 and C-term
Target	G9A/EHMT2
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q96KQ7
NCBI	NP_006700.2
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.