

Product Datasheet

FOG1/ZFPM1 Antibody (orb1518183)

Description	FOG1/ZFPM1 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 175 and 225
Target	FOG1/ZFPM1
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	250 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q8IX07
NCBI	NP_722520.2
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.