

Product Datasheet

Rab5 antibody [Out of stock] (orb151819)

Description

Rabbit polyclonal to Rab5 (ATTO-488).

of small guanosine tr...

Species/Host

Rabbit

Reactivity

Bovine, Human, Monkey, Mouse, Rat

Conjugation

ATTO-488

Tested

ICC, IHC, WB

Applications
Immunogen

Human Rab5 synthetic peptide conjugated to KLH; identical to dog Rab5 sequence over the residues

Target

Rab5

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

1 µl/ml of SPC-168 was sufficient for detection of Rab5 in 15 µg of HeLa lysate by ECL immunoblot analysis using Donkey anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

26kDa

Uniprot ID
[P20339](#)
NCBI
[NP_004153.2](#)
Entrez
[5868](#)
Dilution Range

WB (1:1000), IHC (1:100), ICC/IF (1:80)

Expiration Date

12 months from date of receipt.