



Product Datasheet

14-3-3 Sigma Antibody (orb1518220)


Description	14-3-3 Sigma Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 198 and 248
Target	14-3-3 Sigma
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	50 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P31947
NCBI	NP_006133.1
Hazard Information	Non-Toxic
Dilution Range	IHC - 1:100 - 1:500; IHC-IF - 1:50 - 1:500
Expiration Date	12 months from date of receipt.



Detection of human 14-3-3 Sigma by immunofluorescence.



Detection of human 14-3-3 Sigma by immunohistochemistry (IHC).



Detection of human 14-3-3 Sigma by immunohistochemistry (IHC).