

Product Datasheet

TNF-R1 antibody [Out of stock] (orb151852)

Description

Rabbit polyclonal antibody against TNFR1

Species/Host	Rabbit
Reactivity	Bovine, Canine, Human, Monkey, Mouse, Rabbit, Rat
Conjugation	ATTO-390
Tested Applications	ICC, IHC, IP, WB
Immunogen	Peptide corresponding to AA 20-43 of the mouse TNF-R1 sequence, identical to rat and human over those residues
Target	TNFR1
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SPC-170 was sufficient for detection of TNFR1 in 20 µg of Hela lysate by colorimetric immunoblot analysis using Goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	55kDa
Uniprot ID	P19438
NCBI	P19438
Entrez	8666
Dilution Range	WB (1:1000), IHC (1:100), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.