

Product Datasheet

p38 antibody [Out of stock] (orb151869)

Description

Rabbit polyclonal to p38 (ATTO-390). The

comprises a family ...

Species/Host

Rabbit

Reactivity

Bovine, Canine, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep

Conjugation

ATTO-390

Tested

ICC, IHC, IP, WB

Applications

Immunogen

A 20 residue synthetic peptide based on the human p38 with the cysteine residue added and coupled to KLH

Target

p38

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-172 was sufficient for detection of p38 in 20 µg of HeLa cell lysate by ECL immunoblot analysis.

Clonality

Polyclonal

MW

43kDa

Uniprot ID

[Q16539](#)

NCBI

[NP_001306.1](#)

Entrez

[1432](#)

Dilution Range

WB (1:1000), ICC/IF (1:100), IP (1:250)

Expiration Date

12 months from date of receipt.