

Product Datasheet

Histone H3, Phospho (S10) Antibody (orb1518730)

Description	Histone H3, Phospho (S10) Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ICC, IHC
Immunogen	Surrounding serine 10
Target	Histone H3, Phospho (S10)
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	20 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P68431
NCBI	NP_003520.1
Dilution Range	IHC - 1:100 - 1:500; ICC-IF - 1:100 - 1:500
Expiration Date	12 months from date of receipt.