

Product Datasheet

p38 antibody (orb151880)

Description

Rabbit polyclonal to p38 (FITC). The MAPK

Species/Host

Rabbit

Reactivity

Bovine, Canine, Gallus, Guinea pig, Hamster, Human, Monkey, Mouse, Porcine, Rabbit, Rat, Sheep

Conjugation

FITC

Tested

ICC, IHC, IP, WB

Applications
Immunogen

A 20 residue synthetic peptide based on the human p38 with the cysteine residue added and coupled to KLH

Target

p38

Preservatives

PBS pH 7.4, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

A 1:1000 dilution of SPC-172 was sufficient for detection of p38 in 20 µg of HeLa cell lysate by ECL immunoblot analysis.

Clonality

Polyclonal

MW

43kDa

Uniprot ID
[Q16539](#)
NCBI
[NP_001306.1](#)
Entrez
[1432](#)
Dilution Range

WB (1:1000), ICC/IF (1:100), IP (1:250)

Expiration Date

12 months from date of receipt.