



Product Datasheet

53BP1 Antibody (orb1518840)

Description	53BP1 Antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	IHC
Immunogen	Between 350 and 400
Target	53BP1
Preservatives	Tris-buffered Saline containing 0.1% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	Q12888
NCBI	NP_005648.1
Dilution Range	IHC - 1:100 - 1:500
Expiration Date	12 months from date of receipt.



Detection of human 53BP1 by immunohistoc...



Detection of mouse 53BP1 by immunohistoc...