

Product Datasheet

SARS-CoV-2 IgG ELISA Kit (For Research Use Only) (orb1518853)

Description	SARS-CoV-2 IgG ELISA Kit (For Research Use Only)
Reactivity	Human
Tested Applications	ELISA
Target	SARS-CoV-2 IgG
Storage	2 - 8°C
Note	For research use only
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	Please enquire.