

## Product Datasheet

### Human IgG ELISA Kit (orb1518862)

|                            |                                                                                                                      |
|----------------------------|----------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human IgG ELISA Kit                                                                                                  |
| <b>Reactivity</b>          | Human                                                                                                                |
| <b>Tested Applications</b> | ELISA                                                                                                                |
| <b>Target</b>              | Human IgG                                                                                                            |
| <b>Storage</b>             | 2 - 8°C                                                                                                              |
| <b>Note</b>                | For research use only                                                                                                |
| <b>Dilution Range</b>      | WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF: |
| <b>Expiration Date</b>     | Please enquire.                                                                                                      |