

Product Datasheet

Human Hemoglobin F ELISA Quantitation Set (orb1518865)

Description	Human Hemoglobin F ELISA Quantitation Set
Reactivity	Human
Tested Applications	ELISA
Target	Human Hemoglobin F
Storage	2 - 8°C
Note	For research use only
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP- Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	Please enquire.