

Product Datasheet

Chicken IgA ELISA Kit (orb1518870)

Description	Chicken IgA ELISA Kit
Reactivity	Gallus
Tested Applications	ELISA
Target	Chicken IgA
Storage	2 - 8°C
Note	For research use only
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	Please enquire.