

Product Datasheet

Bovine Lactoferrin ELISA Kit (orb1518875)

Description	Bovine Lactoferrin ELISA Kit
Reactivity	Bovine
Tested Applications	ELISA
Target	Bovine Lactoferrin
Storage	2 - 8°C
Note	For research use only
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	Please enquire.