

Product Datasheet

Pig IgG ELISA Kit (orb1518894)

Description	Pig IgG ELISA Kit
Reactivity	Porcine
Tested Applications	ELISA
Target	Pig IgG
Storage	2 - 8°C
Note	For research use only
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	Please enquire.