

Product Datasheet

Mouse IgG-Heavy and Light Chain Antibody (orb1518899)

Description	Mouse IgG-Heavy and Light Chain Antibody
Species/Host	Goat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	FC, ICC, IF, IHC
Target	Mouse IgG-Heavy and Light Chain
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	IHC - 1:50 - 1:500; ICC - 1:50 - 1:500; IF - 1:50 - 1:500
Expiration Date	12 months from date of receipt.