

Product Datasheet

Mouse IgG-Heavy and Light Chain Antibody (orb1518939)

Description	Mouse IgG-Heavy and Light Chain Antibody
Species/Host	Goat
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Target	Mouse IgG-Heavy and Light Chain
Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Form/Appearance	F(ab') ₂
Concentration	0.5 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: F(ab') ₂
Clonality	Polyclonal
Dilution Range	WB - 1:1,000 - 1:20,000; IHC - 1:200 - 1:2,000; ICC - 1:200 - 1:2,000; ELISA - 1:1,000 - 1:20,000; for coating plates 1:50 - 1:250
Expiration Date	12 months from date of receipt.