

## Product Datasheet

### Mouse IgG2c Antibody (orb1518945)

|                            |                                                                                                                                    |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Mouse IgG2c Antibody                                                                                                               |
| <b>Species/Host</b>        | Goat                                                                                                                               |
| <b>Reactivity</b>          | Mouse                                                                                                                              |
| <b>Conjugation</b>         | Unconjugated                                                                                                                       |
| <b>Tested Applications</b> | ELISA, ICC, IHC, WB                                                                                                                |
| <b>Target</b>              | Mouse IgG2c                                                                                                                        |
| <b>Preservatives</b>       | Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide                                                                      |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                          |
| <b>Concentration</b>       | 0.5 mg/ml                                                                                                                          |
| <b>Storage</b>             | 2 - 8°C                                                                                                                            |
| <b>Note</b>                | For research use only                                                                                                              |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                  |
| <b>Clonality</b>           | Polyclonal                                                                                                                         |
| <b>Dilution Range</b>      | WB - 1:1,000 - 1:20,000; IHC - 1:200 - 1:2,000; ICC - 1:200 - 1:2,000; ELISA - 1:1,000 - 1:20,000; for coating plates 1:50 - 1:250 |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                    |