

Product Datasheet

Mouse Albumin Antibody (orb1518946)

Description	Mouse Albumin Antibody
Species/Host	Goat
Reactivity	Mouse
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC, WB
Target	Mouse Albumin
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.05% Pro-Clean 400
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P07724
Dilution Range	WB - 1:1,000 - 1:30,000; IHC - 1:100 - 1:250; ICC - 1:100 - 1:250; ELISA - 1:10,000 - 1:50,000
Expiration Date	12 months from date of receipt.