

Product Datasheet

Mouse IgA Antibody (orb1519095)

Description	Mouse IgA Antibody
Species/Host	Goat
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC, WB
Target	Mouse IgA-Heavy Chain
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	WB - 1:10,000 - 1:200,000; IHC - 1:250 - 1:2,500; ICC - 1:100 - 1:500; ELISA - 1:10,000 - 1:200,000
Expiration Date	12 months from date of receipt.