

Product Datasheet

ApoB-100 Antibody (orb1519135)

Description	ApoB-100 Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, WB
Immunogen	Human Apolipoprotein B100
Target	ApoB-100
Preservatives	75 mM Sodium Phosphate, 75 mM NaCl, 0.5 mM EDTA, 0.02% NaN ₃ , pH 7.2
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Uniprot ID	P04114
NCBI	NP_000375.2
Dilution Range	WB - 1:1,000 - 1:80,000; ELISA - 1:1,000 - 1:80,000
Expiration Date	12 months from date of receipt.