

## Product Datasheet

### ApoB-100 Antibody (orb1519136)

|                            |                                                                                  |
|----------------------------|----------------------------------------------------------------------------------|
| <b>Description</b>         | ApoB-100 Antibody                                                                |
| <b>Species/Host</b>        | Goat                                                                             |
| <b>Reactivity</b>          | Human                                                                            |
| <b>Conjugation</b>         | Unconjugated                                                                     |
| <b>Tested Applications</b> | ELISA, WB                                                                        |
| <b>Immunogen</b>           | Human Apolipoprotein B100                                                        |
| <b>Target</b>              | ApoB-100                                                                         |
| <b>Preservatives</b>       | 75 mM Sodium Phosphate, 75 mM NaCl, 0.5 mM EDTA, 0.02% NaN <sub>3</sub> , pH 7.2 |
| <b>Form/Appearance</b>     | Whole IgG                                                                        |
| <b>Concentration</b>       | 1 mg/ml                                                                          |
| <b>Storage</b>             | 2 - 8°C                                                                          |
| <b>Note</b>                | For research use only                                                            |
| <b>Application notes</b>   | Format: Whole IgG                                                                |
| <b>Clonality</b>           | Polyclonal                                                                       |
| <b>Uniprot ID</b>          | <b>P04114</b>                                                                    |
| <b>NCBI</b>                | <b>NP_000375.2</b>                                                               |
| <b>Dilution Range</b>      | WB - 1:1,000 - 1:80,000; ELISA - 1:1,000 - 1:80,000                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                  |