


Product Datasheet

CD3E, Cytokeratin, CD8 alpha, CD68, Ki-67, PD-L1 Antibody (orb1519155)

Description	CD3E, Cytokeratin, CD8 alpha, CD68, Ki-67, PD-L1 Antibody
Reactivity	Human
Conjugation	Unconjugated
Target	CD3E, Cytokeratin, CD8 alpha, CD68, Ki-67, PD-L1
Storage	2 - 8°C
Note	For research use only
Clonality	Recombinant
Uniprot ID	P46013, P01732, P08779, P19012, P08727, P34810, P07766, P13645, Q9NZQ7, P02533
NCBI	P02533, P07766, P08727, Q9NZQ7, P13645, P19012, P01732, P46013, P08779, P34810
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	12 months from date of receipt.



Detection of human CD3 (cyan) , CD8 (gre...