

Product Datasheet

CD3E, CD68, CD20 Antibody (orb1519157)

Description	CD3E, CD68, CD20 Antibody
Reactivity	Human
Conjugation	Unconjugated
Target	CD3E, CD68, CD20
Storage	2 - 8°C
Note	For research use only
Clonality	Recombinant
Uniprot ID	P11836, P34810, P07766
NCBI	NP_000724.1, P34810, P11836
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	12 months from date of receipt.