

Product Datasheet

Human IgG-Heavy and Light Chain Antibody (orb1519204)

Description	Human IgG-Heavy and Light Chain Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Biotin
Tested Applications	ELISA, ICC, IHC, WB
Target	Human IgG-Heavy and Light Chain
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	WB - 1:20,000 - 1:400,000; IHC - 1:500 - 1:5,000; ICC - 1:200 - 1:1,000; ELISA - 1:20,000 - 1:400,000
Expiration Date	12 months from date of receipt.