

Product Datasheet

Human IgG-Heavy and Light Chain Antibody (orb1519205)

Description	Human IgG-Heavy and Light Chain Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	AP
Tested Applications	ELISA, ICC, IHC, WB
Target	Human IgG-Heavy and Light Chain
Preservatives	50 mM HEPES pH 7.1, 0.1 M NaCl, 1 mM MgCl ₂ , 0.1 mM ZnCl ₂ containing 0.2% rAlbumin and 0.09% NaN ₃
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Storage	2 - 8°C
Note	For research use only
Clonality	Polyclonal
Dilution Range	WB - 1:2,500 - 1:25,000; IHC - 1:100 - 1:1,000; ICC - 1:100 - 1:1,000; ELISA - 1:2,500 - 1:25,000
Expiration Date	12 months from date of receipt.