

## Product Datasheet

### HSP90 beta antibody [Out of stock] (orb151921)

|                            |                                                                                                                             |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal to Hsp90 beta antibody (ATTO-488)                                                                         |
| <b>Species/Host</b>        | Rabbit                                                                                                                      |
| <b>Conjugation</b>         | ATTO-488                                                                                                                    |
| <b>Tested Applications</b> | ELISA, ICC, IHC, WB                                                                                                         |
| <b>Immunogen</b>           | Full length protein HSP90                                                                                                   |
| <b>Target</b>              | HSP90 beta                                                                                                                  |
| <b>Preservatives</b>       | PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated                                     |
| <b>Concentration</b>       | 1 mg/ml                                                                                                                     |
| <b>Storage</b>             | Conjugated antibodies should be stored at 4°C                                                                               |
| <b>Note</b>                | For research use only                                                                                                       |
| <b>Application notes</b>   | A 1:1000 dilution of SPC-177 was sufficient for detection of HSP90 in 20 µg of HeLa cell lysate by ECL immunoblot analysis. |
| <b>Clonality</b>           | Polyclonal                                                                                                                  |
| <b>MW</b>                  | 90kDa                                                                                                                       |
| <b>Uniprot ID</b>          | <b>P08238</b>                                                                                                               |
| <b>NCBI</b>                | <b>NP_031381.2</b>                                                                                                          |
| <b>Entrez</b>              | <b>3326</b>                                                                                                                 |
| <b>Dilution Range</b>      | WB (1:1000), IHC (1:100), ICC/IF (1:120)                                                                                    |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                             |