

## Product Datasheet

### Human IgM Antibody (orb1519215)

|                            |                                                                                                                                       |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Human IgM Antibody                                                                                                                    |
| <b>Species/Host</b>        | Goat                                                                                                                                  |
| <b>Reactivity</b>          | Human                                                                                                                                 |
| <b>Conjugation</b>         | AP                                                                                                                                    |
| <b>Tested Applications</b> | ELISA, ICC, IHC, WB                                                                                                                   |
| <b>Target</b>              | Human IgM                                                                                                                             |
| <b>Preservatives</b>       | 50 mM HEPES pH 7.1, 0.1 M NaCl, 1 mM MgCl <sub>2</sub> , 0.1 mM ZnCl <sub>2</sub> containing 0.2% rAlbumin and 0.09% NaN <sub>3</sub> |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                             |
| <b>Concentration</b>       | 0.5 mg/ml                                                                                                                             |
| <b>Storage</b>             | 2 - 8°C                                                                                                                               |
| <b>Note</b>                | For research use only                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                            |
| <b>Dilution Range</b>      | WB - 1:2,500 - 1:25,000; IHC - 1:100 - 1:1,000; ICC - 1:100 - 1:1,000; ELISA - 1:2,500 - 1:25,000                                     |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                       |