

Product Datasheet

Human IgM Antibody (orb1519216)

Description	Human IgM Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IHC, WB
Target	Human IgM
Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	0.5 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	WB - 1:1,000 - 1:20,000; IHC - 1:200 - 1:2,000; ICC - 1:200 - 1:2,000; ELISA - 1:1,000 - 1:20,000; for coating plates 1:50 - 1:250
Expiration Date	12 months from date of receipt.