

Product Datasheet

Rabbit anti-Human Alpha-1-Antitrypsin Antibody Affinity Purified (orb1519217)

Description	Rabbit anti-Human Alpha-1-Antitrypsin Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IEP
Target	Human Alpha-1-Antitrypsin
Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	ELISA - 1:1,000 - 1:30,000; for coating plates 1:100 - 1:50; IEP - Neat
Expiration Date	12 months from date of receipt.