

Product Datasheet

Human IgG+IgM+IgA-Heavy and Light Chain Antibody (orb1519223)

Description	Human IgG+IgM+IgA-Heavy and Light Chain Antibody
Species/Host	Goat
Reactivity	Human
Conjugation	HRP
Tested Applications	ELISA, ICC, IHC, WB
Target	Human IgG+IgM+IgA-Heavy and Light Chain
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.05% Pro-Clean 400
Form/Appearance	Whole IgG
Concentration	1 mg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Polyclonal
Dilution Range	WB - 1:5,000 - 1:50,000; IHC - 1:200 - 1:5,000; ICC - 1:200 - 1:5,000; ELISA - 1:10,000 - 1:100,000
Expiration Date	12 months from date of receipt.