

Product Datasheet

HSP70 antibody (orb151947)

Description

Chicken polyclonal to Hsp70 (Biotin). Hsp70

	hsps (hsp70s...
Species/Host	Gallus
Reactivity	Human, Mouse, Rat
Conjugation	Biotin
Tested Applications	ELISA, ICC, WB
Immunogen	Full length Human HSP70
Target	HSP70
Preservatives	PBS pH 7.4, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	1 µg/ml of SMC-178 was sufficient for detection of HSP70 in 20 µg of heat shocked HeLa cell lysate by colorimetric immunoblot analysis using Goat anti-chicken IgY:HRP as the secondary antibody.
Isotype	IgY
Clonality	Polyclonal
MW	70kDa
Uniprot ID	P0DMV8 , P0DMV9
NCBI	NP_005336.3
Entrez	3303
Dilution Range	WB (1:1000), ICC/IF (1:120)
Expiration Date	12 months from date of receipt.