

Product Datasheet

S100A9 Antibody (BSA Free) (orb1519481)

Description	S100A9 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, WB
Immunogen	between 64 and 114
Target	S100A9
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P06702
NCBI	NP_002956.1
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	12 months from date of receipt.