

Product Datasheet

SARS-CoV Spike RBD Antibody (orb1519485)

Description	SARS-CoV Spike RBD Antibody
Species/Host	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	SARS-CoV Spike RBD
Target	SARS-CoV Spike RBD
Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P0DTC2
NCBI	YP_009724390.1
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP-PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	12 months from date of receipt.