

Product Datasheet

GRP78 (Bip) antibody [Out of stock] (orb151954)

Description

Rabbit polyclonal to GRP78 (ATTO-390).

	kDa glucose regulated p...
Species/Host	Rabbit
Reactivity	Canine, Drosophila, Hamster, Human, Mouse, Rat
Conjugation	ATTO-390
Tested Applications	ELISA, ICC, IHC, WB
Immunogen	Full length human GRP78 (Bip) his tagged at the N terminus
Target	GRP78
Preservatives	PBS, 50% glycerol, 0.09% sodium azide *Storage buffer may change when conjugated
Concentration	1 mg/ml
Storage	Conjugated antibodies should be stored at 4°C
Note	For research use only
Application notes	0.5 µg/ml of SPC-180 was sufficient for detection of Grp78 in 10 µg of rat tissue lysate by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.
Clonality	Polyclonal
MW	78kDa
Uniprot ID	P11021
NCBI	NP_005338.1
Entrez	3309
Dilution Range	WB (1:2000), ICC/IF (1:100)
Expiration Date	12 months from date of receipt.