

Product Datasheet

OX40L/TNFSF4 Antibody (orb1519580)

Description	OX40L/TNFSF4 Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	between 1 and 50
Target	OX40L/TNFSF4
Preservatives	Phosphate Buffered Saline (PBS) containing 0.2% rAlbumin and 0.09% Sodium Azide
Form/Appearance	Whole IgG
Concentration	100 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	P23510
NCBI	NP_003317.1
Dilution Range	WB - 1:1,000
Expiration Date	12 months from date of receipt.