

Product Datasheet

Fibroblast Activation Protein alpha/FAP Antibody (BSA Free) (orb1519621)

Description	Fibroblast Activation Protein alpha/FAP Antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, IHC, IP, WB
Immunogen	residues 26-760 (ECD)
Target	Fibroblast Activation Protein alpha/FAP
Preservatives	Borate Buffered Saline (BBS) pH 8.2 with 0.09% Sodium Azide, rAlbumin-Free
Form/Appearance	Whole IgG
Concentration	1000 µg/ml
Storage	2 - 8°C
Note	For research use only
Application notes	Format: Whole IgG
Clonality	Recombinant
Uniprot ID	Q12884
NCBI	NP_004451.2
Dilution Range	WB: , IP: , IHC: , ICC: , IHC-IF: , ICC-IF: , ChIP: , ChIP-chip: , ChIP- PCR: , ChIP-Seq: , ELISA: , IEP: , DD: , IF:
Expiration Date	12 months from date of receipt.