

Product Datasheet

GRP78 (Bip) antibody (orb151963)

Description

Rabbit polyclonal to GRP78 (APC). GRP78 is a

regulated protei...

Species/Host

Rabbit

Reactivity

Canine, Drosophila, Hamster, Human, Mouse, Rat

Conjugation

APC

Tested

ELISA, ICC, IHC, WB

Applications
Immunogen

Full length human GRP78 (Bip) his tagged at the N terminus

Target

GRP78

Preservatives

PBS, 50% glycerol, 0.09% sodium azide
*Storage buffer may change when conjugated

Concentration

1 mg/ml

Storage

Conjugated antibodies should be stored at 4°C

Note

For research use only

Application notes

0.5 µg/ml of SPC-180 was sufficient for detection of Grp78 in 10 µg of rat tissue lysate by ECL immunoblot analysis using goat anti-rabbit IgG:HRP as the secondary antibody.

Clonality

Polyclonal

MW

78kDa

Uniprot ID

[P11021](#)

NCBI

[NP_005338.1](#)

Entrez

[3309](#)

Dilution Range

WB (1:2000), ICC/IF (1:100)

Expiration Date

12 months from date of receipt.