



## Product Datasheet

### UBQLN4/CIP75/Ubiquilin 4 Antibody (orb1519630)

|                            |                                                                                                                                                                                                                                                            |
|----------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | UBQLN4/CIP75/Ubiquilin 4 Antibody                                                                                                                                                                                                                          |
| <b>Species/Host</b>        | Rabbit                                                                                                                                                                                                                                                     |
| <b>Reactivity</b>          | Human                                                                                                                                                                                                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                                                                                                                                                                                                               |
| <b>Tested Applications</b> | FC, ICC, IHC, IP, WB                                                                                                                                                                                                                                       |
| <b>Immunogen</b>           | residues 2-601 (full length)                                                                                                                                                                                                                               |
| <b>Target</b>              | UBQLN4/CIP75/Ubiquilin 4                                                                                                                                                                                                                                   |
| <b>Preservatives</b>       | Borate Buffered Saline (BBS) pH 8.2 with 0.1% rAlbumin and 0.09% Sodium Azide                                                                                                                                                                              |
| <b>Form/Appearance</b>     | Whole IgG                                                                                                                                                                                                                                                  |
| <b>Concentration</b>       | 50 µg/ml                                                                                                                                                                                                                                                   |
| <b>Storage</b>             | 2 - 8°C                                                                                                                                                                                                                                                    |
| <b>Note</b>                | For research use only                                                                                                                                                                                                                                      |
| <b>Application notes</b>   | Format: Whole IgG                                                                                                                                                                                                                                          |
| <b>Clonality</b>           | Recombinant                                                                                                                                                                                                                                                |
| <b>Uniprot ID</b>          | <b>Q9NRR5</b>                                                                                                                                                                                                                                              |
| <b>NCBI</b>                | <b>NP_064516.2</b>                                                                                                                                                                                                                                         |
| <b>Dilution Range</b>      | WB - 1:1,000; IP - 20 µl/1 mg lysate; IHC - 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE tissue sections.; ICC - 1:100 to 1:500. Epitope retrieval with citrate buffer pH6.0 is recommended for FFPE cell sections. |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                                                                                                                                                            |



Detection of human UBQLN4 by western blo...



Detection of human UBQLN4/CIP75/Ubiquilin...